




## POLLED SIRE

	HOMOZYGOUS PP OR POS	HETEROZYGOUS POC OR PO OR P
HOMOZYGOUS PP OR POS	100% POLLED	100% POLLED
HETEROZYGOUS POC OR PO OR P	100% POLLED	75% POLLED
NON- POLLED	100% POLLED	50% POLLED

DAM




## TOP POLLED GENOMICS

Semex has 25 Polled Sires with GTPI

greater than 2700


CDN, AUG\*14



Test through Elevate®

Select with SemexWorks™

Complete your strategy with OptiMate™



While the costs and regulatory burden on dehorning is increasing, so is the merit and

number of polled Holstein and Jersey sires

SPURLOCK, 2014

 **Polled**